

Four-sided Applicator (Frame-Type)

The Frame-Type Four-Sided Applicators are designed for the production of uniform films of paint, adhesives and similar products on plane substrates. They combine the accuracy of fixed applicators with the versatility of multiple gap choice in one unit. These applicators are suitable for use of aqueous, acid, and alkaline products.

It is made of high-grade, corrosion-resistant stainless steel 440C.



Main Technical Parameters:

- ★ Wet Film Width of: 80mm
- ★ The total length: 100mm

★ Ordering information:

BGD 204/1: (25µm, 50µm, 75µm, 100µm)

BGD 204/2: (50µm, 100µm, 150µm, 200µm)

BGD 204/3: (30µm, 60µm, 90µm, 120µm)

BGD 204/4: (5µm, 10µm, 15µm, 20µm)