


# Fume Hood



## TECHNICAL PARAMETERS

Model	FH1000(A)	FH1200(A)	FH1500(A)	FH1800(A)	
External Size (W*D*H)	1040*750*2200mm	1240*800*2200mm	1540*800*2200mm	1840*800*2200mm	
Internal Size (W*D*H)	820*620*770mm	1020*670*770mm	1320*670*770mm	1620*670*770mm	
Work Surface Height	850mm				
Max Opening	450mm				
Air Velocity	0.3-0.8m/s				
Noise	≤60dB				
Exhaust Duct	PVC, Standard length: 4 meters. Φ300mm				
Pipe Strap	1 pc				
Fluorescent Lamp	14W*1	21W*1	28W*1	36W*1	
UV Lamp	Emission of 253.7 nanometers for most efficient decontamination				
Blower	Built-in centrifugal blower; Speed adjustable.				
Front Window	5mm toughened glass; Motorized; Height adjustable.				
Power Supply	110/220V±10%, 60/50Hz				
Consumption	400W	400W	500W	500W	
Material	Exterior	Cold-rolled steel with anti-bacteria powder coating.			
	Interior	High grade melamine board with good acid and alkali resistance function			
	Work Table	Solid chemical resistant physiochemical board			
Standard Accessory	UV lamp*2, Fluorescent lamp, Water tap, Gas tap, Water sink, Base cabinet				
	Total load of 2 waterproof sockets: 1000W 4-meter PVC exhaust duct, pipe strap				
Optional Accessory	Active carbon filter or HEPA filter				
Gross Weight	230kg	270kg	370kg	420kg	
Package	Main Body	1180*900*1750 mm	1380*940*1790mm	1680*940*1790mm	1980*940*1790mm
	Base Cabinet	1180*900*900mm	1380*940*900mm	1680*940*900mm	1980*940*900mm



**LED Display**  
Digital LED display: Airflow Level



**Water&Gas Tap**



**Water&Gas Remote Control**



**Double Layers**  
1.2mm sheet metal surface,  
melamine board inward.