

# Advantage-Lab

The best priced laboratory equipment

## FULLY AUTOCLAVABLE

High Performance Micro Pipette

- Fully Autoclavable
- Light Weight & Ergonomic Design
- Click stop volume setting  
digital system
- Calibration conforms to DIN 12650  
& EN ISO 8655 standards  
ensuring high accuracy & precision



## FULLY AUTOCLAVABLE

### High Performance Micro Pipette

#### *Variable Volume*

Cat No	Volumes	Increments
AV-100	0.5 - 10 µl	0.1 µl
AV-200	2 - 20 µl	0.1 µl
AV-300	5 - 50 µl	0.5 µl
AV-400	10 - 100 µl	0.5 µl
AV-500	20 - 200 µl	1 µl
AV-600	100 - 1000 µl	5 µl
AV-700	500 - 5000 µl	50 µl
AV-800	1 - 10 ml	0.1 ml

#### *Fix Volume*

Cat No	Volume
AF-0005	5 µl
AF-0010	10 µl
AF-0020	20 µl
AF-0025	25 µl
AF-0050	50 µl
AF-0100	100 µl
AF-0200	200 µl
AF-0250	250 µl
AF-0500	500 µl
AF-1000	1000 µl
AF-2000	2000 µl
AF-5000	5000 µl



#### Advantage-Lab GmbH

Otto-Röhm-Strasse 69 - 64293 Darmstadt, Germany

Tel.+49 61 51 82 87 812 - Fax.+49 61 51 82 87 859

Web : [www.advantage.com](http://www.advantage.com)

E-mail : [info@advantage-lab.com](mailto:info@advantage-lab.com)

