

Our LR range offers high temperature uniformity and stability and is produced with attention to detail as well as user-friendly features. Products are available with different capacities and offer you the best in technology and reliability.

Features on LR range

- ARCTIKO Controller G214
- White or stainless steel exterior
- Stainless steel interior
- Automatic door closing
- Air forced cooling
- Adjustable shelves
- Castors
- Key lock on door
- Reversible door



Available with stainless steel exterior

LR 300



Arctiko G214 Controller
Read page 4-5



Stainless steel wire shelf (optional)



Perforated stainless steel shelf (optional)



Stainless steel drawer (optional)

Our universal shelf and drawer guides allow you to mix shelves and drawers for flexible storage in each of our biomedical refrigerators and freezers.



OPTIONAL SPARK FREE INTERIOR
ATEX approved according to EN 60079-15
category 3, zone 2. CE EX 11-3 G



LR 1400



LR 660-2



LR 100

Under counter model

Model	LR 100	LR 270-2	LR 300	LR 500	LR 660-2	LR 700	LR 900	LR 1400
Item code (white exterior)	DAI 0280	DAI 0230	DAI 0220	DAI 0240	DAI 0260	DAI 0250	DAI 0335	DAI 0270
Item code (stainless steel exterior)	DAI 0281	DAI 0231	DAI 0221	DAI 0241	DAI 0261	DAI 0251	DAI 0336	DAI 0271
Item code (spark free interior – ATEX)	DAI 0280/AT	DAI 0230/AT	DAI 0220/AT	DAI 0240/AT	DAI 0260/AT	DAI 0250/AT	DAI 0335/AT	DAI 0270/AT
Dimensions exterior WxDxH (mm)	610x659x790	520x700x1997	520x700x1997	620x860x1997	720x860x1997	720x860x1997	1021x860x1997	1440x860x1997
Capacity (L)	94	161/161	346	515	288/288	618	866	1361
Temperature range (°C)	+1/+10	+1/+10	+1/+10	+1/+10	+1/+10	+1/+10	+1/+10	+1/+10
Power supply (V)	230	161 / 161	230	230	230	230	230	230
Frequency (Hz)	50/60	50/60	50/60	50/60	50/60	50/60	50/60	50/60

We reserve the right to change the specifications without notice. Subject to confirmation, availability and errors.

Offering outstanding cooling performance, glass door freezers from Arctiko allow the best overview of the inside of the units. The advanced Arctiko controller ensures reliable temperatures and comes equipped with many alarms as standard.

Features on PR range

- ARCTIKO Controller G214
- White or stainless steel exterior
- Stainless steel interior
- Automatic door closing
- Air forced cooling
- Adjustable shelves
- Castors
- Key lock on door
- Reversible door
- Glass door
- LED light



PR 500

Available with stainless steel exterior



Arctiko G214 Controller
Read page 4-5



Stainless steel wire shelf (optional)



Perforated stainless steel shelf (optional)



Stainless steel drawer (optional)

Our universal shelf and drawer guides allow you to mix shelves and drawers for flexible storage in each of our biomedical refrigerators and freezers.



OPTIONAL SPARK FREE INTERIOR
ATEX approved according to EN 60079-15
category 3, zone 2. CE EX 11-/3 G



PR 1400



PR 700



PR 100

Under counter model

Model	PR 100	PR 300	PR 500	PR 700	PR 900	PR 1400
Item code (white exterior)	DAI 0895	DAI 0885	DAI 0850	DAI 0860	DAI 0810	DAI 0870
Item code (stainless steel exterior)	DAI 0896	DAI 0886	DAI 0851	DAI 0861	DAI 0811	DAI 0871
Item code (spark free interior - ATEX)	DAI 0895/AT	DAI 0885/AT	DAI 0850/AT	DAI 0860/AT	DAI 0810/AT	DAI 0870/AT
Dimensions exterior WxDxH (mm)	610x659x790	520x690x1997	620x860x1997	720x860x1997	1020x860x1997	1440x860x1997
Capacity (L)	94	352	523	628	878	1381
Temperature range (°C)	+1/+10	+1/+10	+1/+10	+1/+10	+1/+10	+1/+10
Power supply (V)	230	230	230	230	230	230
Frequency (Hz)	50/60	50/60	50/60	50/60	50/60	50/60